
System Requirements - Medical Transcriptionist with HTML Editor

If you are not an “administrator” on your computer system you will need to have registry access given to you. Please contact your Information Services department for access.

Hardware

Minimum 450 MHz Processor - Pentium 3 required, 700+ MHz preferred

- Minimum 128 MB of RAM (256 MB preferred)
- 2 GB of available hard disk space
- 1024 x 768 resolution
- Sound Card (NOTE: 32 bit Sound card with 15 pin Game Port needed if using 15 pin game port-based foot pedal)
- 56k minimum; DSL, Cable, T1 provide improved performance
- Headphones (Headphones plugged into external amplified speakers preferred)
- Keyboard, Mouse, Monitor (17” minimum)
- Foot pedal (see table below)

Operating System	Windows 98	Windows 2000	Windows XP
Foot Pedal/Port	15 Pin connector/ Game Port	15 Pin connector/ Game Port 9 pin connector/ Serial Port	15 Pin connector/ Game Port 9 pin connector/ Serial Port

Software

- IE 6.0 or above is required (see “Updating Configuration Settings for Internet Explorer” on page 4 for browser settings)
- Windows 98, Windows 2000, Windows XP (Service Pack 1)
- Windows Media Player up to version 8
- Adobe Acrobat Reader[®] (For access to documentation)

For free download and use of Adobe Acrobat Reader[®] visit the Adobe web site: <http://www.adobe.com/products/acrobat>

Service

- ISP Service NOTE: Free ISPs are not supported.

System Requirements - Medical Transcriptionist with MS Word Editor

If you are not an “administrator” on your computer system you will need to have registry access given to you. Please contact your Information Services department for access.

Hardware

Minimum 450 MHz Processor - Pentium 3 required, 700+ MHz preferred

- Minimum 256 MB of RAM (512 MB preferred)
- 2 GB of available hard disk space
- 1024 x 768 resolution
- Sound Card (NOTE: 32 bit Sound card with 15 pin Game Port needed if using 15 pin game port-based foot pedal)
- 56k minimum; DSL, Cable, T1 provide improved performance
- Headphones (Headphones plugged into external amplified speakers preferred)
- Keyboard, Mouse, Monitor (17” minimum)

- Foot pedal (see table below)

Operating System	Windows 98	Windows 2000	Windows XP
Foot Pedal/ Port	15 Pin connector/ Game Port	15 Pin connector/Game Port 9 pin connector/ Serial Port	15 Pin connector/ Game Port 9 pin connector/ Serial Port

Software

- IE 6.0 or above is required (see“Updating Configuration Settings for Internet Explorer” on page 4 for browser settings)
- Windows 98, Windows 2000, Windows XP (Service Pack 1)
- Windows Media Player up to version 8
- **Microsoft Word 2002 Required**
- Adobe Acrobat Reader[®] (For access to documentation)

For free download and use of Adobe Acrobat Reader[®] visit the Adobe web site: <http://www.adobe.com/products/acrobat>

Service

- ISP Service NOTE: Free ISPs are not supported.

Updating Configuration Settings for Internet Explorer

The following settings must be configured properly in order for Vianeta products to work correctly:

1. Select the menu item **Tools > Internet Options...**
2. Click **Temporary Internet Files > Settings**.
3. Make sure **Every visit to the page** is checked on. Click **OK**.
4. Select **Privacy > Advanced**.
5. Make sure **Override automatic cookie handling** and **Always allow session cookies** are checked on. Click **OK**.
6. Select **Security > Custom Level**.
7. Scroll down to the Miscellaneous section. Make sure **Display Mixed Content** is set to **Enable**. Click **OK**.
8. Under **Select a Web content zone ...**, click the **Trusted Sites** button. Click the **Sites** button.
9. In the **Add this Web site to the zone:** field, enter the value that is contained in **Server URL** field in the **Tools > Edit Preferences...** dialog. Click **Add**.

Click **OK** to close Internet Options.